coLigne

interbody fusion cages

Package insert
Packungsbeiblatt
Encart
Indicazioni d'uso

coligne

CE 0297

coligne.com

coligne AG

Utoquai 43, CH-8008 Zurich, Switzerland
Phone +41 (0)43 343 8000, Fax +41 (0)43 343 8009

© 2020 coLigne AG. All rights reserved. coLigne, the coLigne symbol and ostaPek are registered trademarks of coLigne AG.

Attention, see instruction for use
Attention, consulter le mode d'emploi
Achtung, siehe Gebrauchsanweisung
Attenzione, fare riferimento alle istruzioni per l'uso

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Interbody fusion cages

Caution

Legal requirements restrict these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required.

Important note to operating surgeon

The coLigne ostaPek® interbody fusion cages are designed to assist in providing an adequate biomechanical environment for fusion and should be used with bone graft. Without bone graft, their use may not be successful, as the cage serves as a mechanical strut for a temporary period while the surrounding bone heals and carries the spinal load. Interbody fusion surgical procedure should only be undertaken after the surgeon has had hands on training in this method, and has become thoroughly knowledgeable about spinal anatomy and biomechanics.

Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implant and possible adverse effects of the surgery. The patient should be instructed to limit postoperative activity as this will reduce the risk of bending, breaking or loosening the implant. The patient must be made aware that implant may bend, break or loosen even though restrictions in activity are followed.

Postoperative evaluation of the fusion and implant status is necessary. After the fusion, implant has not to be removed.

Description

The coLigne ostaPek® interbody fusion cages consist of one and/or more ostaPek® cage in a variety of sizes and is available for anterior cervical approach, anterior, lateral and posterior thoracolumbar approach. Instruments required to implant the device are also available. The system components should not be considered for removal following fusion.

Implant materials

The following materials are available as implant materials in the coLigne ostaPek® interbody fusion cages:

1. ostaPek® (Long Carbon Fiber Reinforced Polymer, LCFRP)

Improper combinations of implants and/or materials can have detrimental effect on the mechanical function and longevity of the implant. Further information on the chemical composition and the mechanical properties of the materials is contained in the corresponding material catalogue sheet or in the corresponding product information, which may be obtained from the local representatives.

Storage and handling

Implants are extremely sensitive to damages. Every small scratch or mark of impact on the surfaces can cause excessive wear and can give rise to complications. Extremely careful handling is strongly recommended. Implants must be stored in the original packaging, unopened. If a loaner or consignment set is used it should be checked for completeness. All items should be checked to ensure functionality and maintenance status prior to use. Damaged, expired or non functioning components should not be used and must be returned to the responsible local representative for repair or replacement. The instruments are subject to a certain degree of wear and have to be considered non-durable materials. Before use, they must be checked for correct functioning, and, if necessary, they must be returned to the responsible local representative for repair or replacement. For additional information refer to the coLigne guidelines for maintenance, cleaning and sterilization of implants and instruments.

Cleaning and decontamination

Instruments and implants must be cleaned and sterilized before being used in the surgical field. In addition, instruments and implants previously used in the surgical theater must be decontaminated, cleaned and sterilized before reuse in the surgical field. Decontamination and cleaning should be performed and documented according to established hospital methods (coLigne recommendation). Instruments might require dismantling prior to decontamination and cleaning. Handling of instruments and implants should be performed with care. Improper use or handling may cause damage and lost of functionality.

Sterilization

The coLigne ostaPek® interbody fusion cages are supplied clean and non sterile. All implants and instruments should be cleaned and sterilized prior to surgery. coLigne recommended practices for in hospital cleaning and sterilization should be followed.

Sterilization testing of components has shown the following recommendations for sterilization are effective:

Method: Steam

Cycle: Pre-vacuum (Wrapped)

Temperature: 134° C/273° F

Exposure Time: 18 minutes

These parameters are validated to sterilize only these devices. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The autoclave must be properly installed, maintained and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

Precautions

1. Implant selection. The coLigne ostaPek® Interbody fusion cages are available in a variety of sizes to insure proper sizing of implanted components. The potential for the success of the fusion is increased by selecting the correct size and shape of the implant. Undersizing of implants can lead to premature failure of the implant.

2. Delayed union or nonunion. The coLigne ostaPek® interbody fusion cages are designed to assist in providing an adequate biomechanical environment for fusion. If a delayed union or nonunion occurs the implant may fail due to fatigue. Patients should be fully informed of the risk of implant failure and the possible consequences.

3. Patient selection. Proper patient selection is critical to the success of the procedure. Only patients who satisfy the criteria set forth under the indications section of this document and who do not have any of the conditions set forth under the contraindications section of this document should be considered for interbody fusion cage surgery using the coLigne ostaPek® interbody fusion cages.

The benefit of spinal fusions using any interbody fusion cage system has not been adequately established in patients with stable spines. In addition, patients who smoke have been shown to have an increased incidence of pseudarthrosis. Patients should be informed about the possible consequences.

4. Single use only. Implants are provided as single use only implants and are not to be reused or re-implanted after first implantation in a patient regardless of an apparent undamaged condition.

5. Handling. Implants should be handled and stored appropriately to protect them from unintentional damage. The surgeon should avoid introducing notches or scratches into the implant surface as these may induce premature failure of the implant. Care must be taken when placing an implant to avoid damage.

6. Patient education. Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implant and possible adverse effects of the surgery. The patient should be instructed to limit postoperative activity, as this will reduce the risk of bent, broken or loose implant. The patient must be made aware that implant may bend, break or loosen even though restriction in activity are followed.

Indications

The coLigne ostaPek® interbody fusion cages are intended for patients with disc disease and are indicated for a spinal fusion of the segment involved.

Contraindications

Use of the coLigne ostaPek® interbody fusion cages is contraindicated when there was recently or there is an active infection, whether systemic or localized to the site of the proposed implantation. Any condition that precludes the possibility of fusion is relative contraindication. This includes but is not limited to: cancer, fever, mental illness, alcoholism or drug abuse, osteoporosis or osteopenia, obesity, pregnancy, foreign body sensitivity, extreme instability.

See also the warnings, precautions and possible adverse affects sections of this insert. The following are specific warnings, precautions and adverse effects which should be understood by the surgeon and explained to the patient. General surgical risks should be explained to the patient prior to surgery.

Warnings

Potential risks identified with the use of the coLigne ostaPek® interbody fusion cages, which may require additional surgery, include:

1. Implant fracture,
2. Loss of fixation,
3. Non-union,
4. Fracture of the vertebra,
5. Neurologic injury,
6. Vascular or visceral injury.

See also precautions and possible adverse events sections of the package insert for a complete list of potential risks.

Possible adverse events

Potential risks identified with the use of the coLigne ostaPek® interbody fusion cages, which may require additional surgery, include but are not limited to:

1. Bending, fracture or loosening of implant,
2. Nonunion or delayed union,
3. Fracture of the vertebra,
4. Neurological, vascular or visceral injury,
5. Allergic reaction to a foreign body,
6. Infection,
7. Decrease in bone density due to stress shielding,
8. Pain, discomfort, or abnormal sensations due to the presence of the device,
9. Nerve damage due to surgical trauma,
10. Subsidence,
11. Bursitis,
12. Dural leak,
13. Paralysis,
14. Death.

Note

Additional surgery may be necessary to correct some of these anticipated adverse reactions.

Limited warranty

coLigne AG products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

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coLigne AG, Utoquai 43, CH-8008 Zurich