coLigne

pedicle fixation systems

Package insert
Packungsbeiblatt
Encart
Indicazioni d'uso

coligne

CE 0297

coligne.com

coligne AG
Utoquai 43, CH-8008 Zurich, Switzerland
Phone +41 (0)43 343 8000, Fax +41 (0)43 343 8009

© 2020 coLigne AG. All rights reserved. coLigne, the coLigne symbol, and ostaPek are registered trademarks of coLigne AG.

Attention, see instruction for use
Attention, consulter le mode d'emploi
Achtung, siehe Gebrauchsanweisung
Attenzione, fare riferimento alle istruzioni per l'uso

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pedicle fixation systems

Caution

Legal requirements restrict these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required.

Important note to operating surgeon

The coLigne pedicle fixation systems are designed to aid in the surgical correction and stabilization of the spine during the development of a solid fusion mass. The systems are intended to be used only in the thoracic, lumbar and sacral levels of the spine. Pedicle fixation should only be undertaken after the surgeon has had hands on training in this method of pedicle fixation and has become thoroughly knowledgeable about spinal anatomy and biomechanics.

Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implants and possible adverse effects of the surgery. The patient should be instructed to limit postoperative activity, as this will reduce the risk of bending, breaking or loosening the implants. The patient must be made aware that implants may bend, break or loosen even though restrictions in activity are followed. Postoperative evaluation of the fusion and implant status is necessary. The surgeon should consider removing the implant once a solid fusion is obtained. The patient must be informed of the potential of this secondary surgical procedure and the associated risks.

Description

The coLigne pedicle fixation systems consist of a variety of shapes and sizes of plates, connectors, screws, optional washers, hooks, nuts, rods and transverse connectors. Rods and plates are made of ostaPek® or titanium alloy and are available in a variety of lengths and contours to approximate individual patients needs. Special self-breaking nuts are provided for plate-rod/connector/screw fixation. Always tighten until the upper part of the nut breaks off to leave the assembly at optimum fixation security. The head part of self-breaking nuts and/or screws should not remain in the patient. After the upper part of the self-breaking nuts and/or screws has been sheared off, further re-tightening is not advised. The fixation nut must be removed and totally replaced if any subsequent re-adjustment is necessary.

The coLigne pedicle fixation systems contain components that are intended to be used only with coLigne system parts. They should not be mixed with any other instrumentation, since other fixation systems may not be compatible nor guaranteed.

Device specific instruments should be used during implantation. Instruments required to retrieve the device are also available. The coLigne pedicle fixation systems implants are temporary devices and should be considered for removal following fusion.

Implant materials

The following materials are available as implant materials in the coLigne pedicle fixation systems:

1. Titanium alloy (ISO 5832-3)
2. ostaPek® (Long Carbon Fiber Reinforced Polymer, LCFRP)

Improper combinations of implants and/or materials can have detrimental effect on the mechanical function and longevity of the implant. Further information on the chemical composition and the mechanical properties of the materials is contained in the corresponding material catalogue sheet or in the corresponding product information, which may be obtained from the local representatives.

Storage and handling

Implants are extremely sensitive to damages. Every small scratch or mark of impact on the surfaces can cause excessive wear and can give rise to complications. Extremely careful handling is strongly recommended. Implants must be stored in the original packaging, unopened. If a loaner or consignment set is used it should be checked for completeness. All items should be checked to ensure functionality and maintenance status prior to use. Damaged, expired or non functioning components should not be used and must be returned to the responsible local representative for repair or replacement. The instruments are subject to a certain degree of wear and have to be considered non-durable materials. Before use, they must be checked for correct functioning, and, if necessary, they must be returned to the responsible local representative for repair or replacement. For additional information refer to the coLigne guidelines for maintenance, cleaning and sterilization of implants and instruments.

Cleaning and decontamination

Instruments and implants must be cleaned and sterilized before being used in the surgical field. In addition, instruments and implants previously used in the surgical theater must be decontaminated, cleaned and sterilized before reuse in the surgical field. Decontamination and cleaning should be performed and documented according to established hospital methods (coLigne recommendation). Instruments might require dismantling prior to decontamination and cleaning. Handling of instruments and implants should be performed with care. Improper use or handling may cause damage and loss of functionality.

Sterilization

The coLigne pedicle fixation system components are supplied clean and non sterile. All implants and instruments should be cleaned and sterilized prior to surgery. coLigne recommended practices for in hospital cleaning and sterilization should be followed. Sterilization testing of components has shown the following recommendations for sterilization are effective:
Method: Steam
Cycle: Pre-vacuum (Wrapped)
Temperature: 134° C/273° F
Exposure Time: 18 minutes

These parameters are validated to sterilize only these devices. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The autoclave must be properly installed, maintained and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

Precautions

1. Implant selection. The coLigne pedicle fixation systems are available in a variety of sizes to insure proper sizing of implanted components. The potential for the success of the fusion is increased by selecting the correct size and shape of the implant. Undersizing of implants can lead to premature failure of the implants. These devices are not intended to be used as the sole support for the spine.

2. Delayed union or nonunion. The coLigne pedicle fixation systems are designed to assist in providing an adequate biomechanical environment for fusion. They are not intended to and must not be used as the sole support for the spine. If a delayed union or nonunion occurs the implant may fail due to fatigue. Patients should be fully informed of the risk of implant failure and the possible consequences.

3. Improper mixing of materials. The coLigne pedicle fixation implants are available in titanium alloy and ostaPek®. It is imperative that titanium alloy does not come into contact in vivo with other dissimilar metals. Accelerated corrosion may occur when two dissimilar metals are in contact within the body environment.

4. Patient selection. Proper patient selection is critical to the success of the procedure. Only patients who satisfy the criteria set forth under the indications section of this document and who do not have any of the conditions set forth under the contraindications section of this document should be considered for pedicle fixation surgery using the coLigne pedicle fixation systems.

The benefit of spinal fusion utilizing any pedicle screw fixation system has not been adequately established in patients with stable spines. In addition, patients who smoke have been shown to have an increased incidence of pseudarthrosis. Patients should be advised and informed about the possible consequences.

5. Single use only. Implants are provided as single use only implants and are not to be reused or reimplanted after first implantation in a patient regardless of an apparent undamaged condition.

6. Handling. Implants should be handled and stored appropriately to protect them from unintentional damage. The surgeon should avoid introducing notches or scratches into the implant surface as these may induce premature failure of the implant. Care must be taken when placing an implant to avoid damage.

7. Patient education. Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implant and possible adverse effects of the surgery. The patient should be instructed to limit postoperative activity, as this will reduce the risk of bent, broken or loose implants. The patient must be made aware that implants may bend break or loosen even though restrictions in activity are followed.

Indications

When used as a pedicle screw fixation system in the non-cervical posterior spine in skeletally mature patients, the coLigne pedicle fixation systems are indicated for one or more of the following indications:

1. Degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies),
2. Degenerative spondylolisthesis with objective evidence of neurologic impairment,
3. Trauma,
4. Dislocation,
5. Scoliosis,
6. Kyphosis,
7. Spinal tumor,
8. Failed previous fusion (pseudarthrosis).

In addition, when used as a pedicle screw fixation system, the coLigne pedicle fixation systems are indicated for skeletally mature patients:

1. Having severe spondylolisthesis (Grades 3 and 4) of the fifth lumbar-first sacral (L5-S1) joint,
2. Who are receiving fusion using bone graft,
3. Who are having the device fixed or attached to the lumbar and sacral spine (L3 and below),
4. Who are having the device removed after the development of a solid fusion mass.

In addition, when used as an anterolateral thoraco-lumbar system, the coLigne spinal fixation systems are indicated for one or more of the following indication:

1. Degenerative disc disease,
2. Trauma,
3. Spinal deformity,
4. Spinal tumor,
5. Spinal stenosis,
6. Degenerative spondylolisthesis,
7. Failed previous fusion (pseudarthrosis).

The coLigne pedicle fixation systems are also a connector and sacral/iliac screw fixation of the non-cervical spine indicated for:

1. Degenerative disc disease,
2. Degenerative spondylolisthesis with objective evidence of neurologic impairment,
3. Trauma,
4. Spinal stenosis,
5. Scoliosis,

6. Kyphosis,
7. Spinal tumor,
8. Failed previous fusion (pseudarthrosis).

Contraindications

The coLigne pedicle fixation systems are not intended for cervical spine use.

Other contraindications include, but are not limited to:

1. Active infections or high risk of infection, local or systemic,
2. Signs of local inflammation,
3. Fever or leukocytosis,
4. Obesity. A patient that is obese or overweight can place loads upon the implants that can cause implant failure, failure of fixation and prevent the formation of a solid bone mass,
5. Pregnancy,
6. Mental illness,
7. Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of severe congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of the white blood count (WBC), or a marked left shift in the WBC differential count,
8. Rapid joint disease, bone absorption, osteopenia, and/or osteoporosis. Osteoporosis and osteopenia are a relative contraindication: they may limit the degree of obtainable correction and the amount of mechanical fixation,
9. Suspected or documented metal allergy or intolerance,
10. Any case not needing a bone graft and fusion,
11. For pedicle screw cases, missing or deformed lumbar pedicles,
12. Any case in which a wrong selection of the implants size to use is performed,
13. Any case in which the implants are misassembled (mixing implants of two different metals as an example),
14. Any case of mixing implants belonging to systems not cleared by coLigne,
15. Any patient having inadequate tissue coverage over the operative site,
16. Any case not described in the indications,
17. Any patient unwilling to cooperate with the postoperative instructions,
18. Use of this type of surgical implant surgery in children or pediatric patients presents particular risks because of bone growth or physical movement specific to this kind of patients. Subsequent re-intervention may be required,
19. Infants with a known hereditary or acquired bone friability or calcification problem, or those with a very short life expectancy, should not be considered for this type of surgery,
20. Substance abuse or senility that precludes the patient from following postoperative precautions to prevent implant failure,
21. Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance,

Contraindications of these systems are consistent with those of other spinal rod and plate systems for fusion. These spinal implant systems are not designed, intended, or sold for uses other than those indicated.

See also the warnings, precautions and possible adverse affects sections of this insert.

The following are specific warnings, precautions and adverse effects, which should be understood by the surgeon and explained to the patient. These warnings are specific to internal fixation devices. General surgical risks should be explained to the patient prior to surgery.

Warnings

Potential risks identified with the use of the coLigne pedicle fixation systems, which may require additional surgery, include:

1. Implant fracture,
2. Metallosis,
3. Loss of fixation,
4. Non-union,
5. Fracture of the vertebra,
6. Neurologic injury,
7. Vascular or visceral injury.

See also precautions and possible adverse events sections of the package insert for a complete list of potential risks.

Possible adverse events

Potential risks identified with the use of the coLigne pedicle fixation systems, which may require additional surgery, include but are not limited to:

1. Early or late bending, disassembly, fracture or loosening of any or all implant(s),
2. Nonunion, delayed union (or pseudarthrosis),
3. Implant migration,
4. Disassembly, bending, loosening, slippage, and/or breakage of instruments peroperatively,
5. Scar formation leading to neurological compromise, nerve compression, or pain,
6. Decrease in bone density, resorption, fracture of the vertebra and/or of the bone graft, at, above or below the treated level, due to, for example, stress shielding. Retro-pulsed graft,
7. Hernia or any disc diseases at, above or below the treated level,
8. Metal sensitivity, allergic reaction and/or foreign body reaction to the implants including possible tumor formation, autoimmune disease, metallosis, and/or scarring,
9. Infection,

- Pressure on the skin possibly resulting in skin breakdown from component parts where there is inadequate tissue coverage over the implant. Implant or graft extrusion through the skin. Wound complications,
- Loss of proper spinal curvature, correction, height, and/or reduction,
- Loss of neurological function, appearance of radiculopathy, dural tears, and/or development of pain. Neurovascular deficits (temporary or permanent) including paralysis or other types of serious injury. Cerebral fluid leakage, meningitis,
- Gastrointestinal, urological, and/or reproductive system compromise, including sterility. Impotency,
- Hemorrhage of blood vessels and/or hematomae,
- Cessation of growth of the fused portion of the spine,
- Discitis, arachnoiditis, and/or other types of inflammation,
- Deep venous thrombosis, thrombophlebitis, and/or pulmonary embolus. Respiratory complications,
- Bone graft donor site complications,
- Inability to resume activities of normal daily living. Loss or increase of mobility or function,
- Change in psychological status,
- Death.

Note


Additional surgery may be necessary to correct some of these anticipated adverse reactions.

Limited warranty

coLigne AG products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

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coLigne AG, Utoquai 43, CH-8008 Zurigo