coLigne

VBR cages

Package insert
Packungsbeiblatt
Encart
Indicazioni d'uso
Insertar paquete

coligne 

CE 0297

coligne AG
Utoquai 43, CH-8008 Zurich, Switzerland
Phone +41 (0)43 343 8000, Fax +41 (0)43 343 8009
www.coligne.com

© 2023 coLigne AG. All rights reserved. coLigne, the coLigne symbol and ostaPek are registered trademarks of coLigne AG.

Attention, see instruction for use
Attention, consulter le mode d'emploi
Achtung, siehe Gebrauchsanweisung
Attenzione, fare riferimento alle istruzioni per l'uso
Atención, ver instrucciones

Single use only, do not reuse
Usage unique, ne pas réutiliser
Zum einmaligen Gebrauch, nicht wiederverwenden
Esclusivamente monouso, non riutilizzabile
Un solo uso, no reutilizar

Reference number
Numéro de référence
Referenz Nummer
Numero di referencia
Número de referencia
Lot number
Numéro de lot
Lot Nummer
Numero di lotto
Número de lote

Non-Sterile
Non-Stérile
Non-Steril
Non-Sterile
No estéril
Manufacturer
Fabriquant
Hersteller
Fabbricante
Fabricante

EU authorized representative
Mandataire UE
EU Bevollmächtigte
Mandatario UE
Representante autorizado EU

English

VBR cages

Caution
Legal requirements restrict these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required.

Important note to operating surgeon
The coLigne ostaPek® vertebral body replacement (VBR) cages are intended for use in the thoracolumbar or cervical spine for partial or total replacement of a collapsed, damaged or unstable vertebral body due to tumor or trauma (i.e., fracture). The coLigne VBR system is compatible with the coLigne pedicle fixation systems in the thoracolumbar spine. In case of

one or multilevel vertebral body reconstruction, posterior fixation of, at least, the two adjacent segments above and below is recommended. The coLigne ostaPek® VBR cage is designed to be used with bone graft. VBR surgical procedure should only be undertaken after the surgeon has had hands on training in this method, and has become thoroughly knowledgeable about spinal anatomy and biomechanics. The coLigne ostapek® VBR cages shall not be impacted during surgery. Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implant and possible adverse effects of the surgery. The patient should be instructed to limit postoperative activity as this will reduce the risk of bending, breaking or loosening the implants. The patient must be made aware that implants may bend, break or loosen even though restrictions in activity are followed. Postoperative evaluation of the fusion and implant status is necessary. After the fusion, implant has not to be removed.

Description
The coLigne VBR cages consist of one and/or more ostaPek® cage in a variety of sizes. Artificial pedicle connecting elements are also available. Instruments required to implant the device are available. The system components should not be considered for removal following fusion.

Implant materials
The following materials are available as implant materials in the coLigne VBR:
1. ostaPek® (Long Carbon Fiber Reinforced Polymer, LCFRP)
2. Titanium alloy (ISO 5832-3)
Improper combinations of implants and/or materials can have detrimental effect on the mechanical function and longevity of the implant. Further information on the chemical composition and the mechanical properties of the materials is contained in the corresponding material catalogue sheet or in the corresponding product information, which may be obtained from the local representatives.

Storage and handling
Implants are extremely sensitive to damages. Every small scratch or mark of impact on the surfaces can cause excessive wear and can give rise to complications. Extremely careful handling is strongly recommended. Implants must be stored in the original packaging, unopened. If a loaner or consignment set is used it should be checked for completeness. All items should be checked to ensure functionality and maintenance status prior to use. Damaged, expired or non functioning implants should not be used and must be returned to the responsible local representative for repair or replacement. The instruments are subject to a certain degree of wear and have to be considered non-durable materials. Before use, they must be checked for correct functioning, and, if necessary, they must be returned to the responsible local representative for repair or replacement. For additional information refer to the coLigne guidelines for maintenance, cleaning and sterilization of implants and instruments.

Cleaning and decontamination
Instruments and implants must be cleaned and sterilized before being used in the surgical field. In addition, instruments and implants previously used in the surgical theater must be decontaminated, cleaned and sterilized before reuse in the surgical field. Decontamination and cleaning should be performed and documented according to established hospital methods (Coligne recommendation). Instruments might require dismantling prior to decontamination and cleaning. Handling of instruments and implants should be performed with care. Improper use or handling may cause damage and lost of functionality.

Sterilization
The coLigne ostaPek® VBR cages are supplied clean and non sterile. All implants and instruments should be cleaned and sterilized prior to surgery. coLigne recommended practices for in hospital cleaning and sterilization should be followed. Sterilization testing of components has shown the following recommendations for sterilization are effective:
Method: Steam
Cycle: Pre-vacuum (Wrapped)
Temperature: 134° C/273° F
Exposure Time: 18 minutes
These parameters are validated to sterilize only these devices. If other products are added to the sterilizer, the recommended parameters are not valid and new cycle parameters must be established by the user. The autoclave must be properly installed, maintained and calibrated. Ongoing testing must be performed to confirm inactivation of all forms of viable microorganisms.

Precautions
1. Implant selection. The coLigne ostaPek® VBR cages are available in a variety of sizes to insure proper sizing of implanted components. The potential for the success of the fusion is increased by selecting the correct size and shape of the implant. Undersizing of implants can lead to premature failure of the implants.

2. Delayed union or nonunion. The coLigne ostaPek® VBR cages are designed to assist in providing an adequate biomechanical environment for fusion. If a delayed union or nonunion occurs the implant may fail due to fatigue. Patients should be fully informed of the risk of implant failure and the possible consequences.
3. Patient selection. Proper patient selection is critical to the success of the procedure. Only patients who satisfy the criteria set forth under the indications section of this document and who do not have any of the conditions set forth under the contraindications section of this document should be considered for VBR cages surgery using the coLigne ostaPek® VBR cages. Patients who smoke have been shown to have an increased incidence of pseudarthrosis. Patients should be informed about the possible consequences.
4. Single use only. Implants are provided as single use only implants and are not to be reused or re-implanted after first implantation in a patient regardless of an apparent undamaged condition.
5. Handling. Implant should be handled and stored appropriately to protect them from unintentional damage. The surgeon should avoid introducing notches or scratches into the implant surface as these may induce premature failure of the implant. Care must be taken when placing an implant to avoid damage.
6. Patient education. Preoperative instructions to the patient are essential. The patient should be made aware of the limitations of the implant and possible adverse effects of the surgery. The patient should be instructed to limit postoperative activity, as this will reduce the risk of bent, broken or loose implants. The patient must be made aware that implants may bend, break or loosen even though restriction in activity are followed.

Indications
The coLigne ostaPek® vertebral body replacement (VBR) cages are intended for use in the thoracolumbar or cervical spine for partial or total replacement of a collapsed, damaged or unstable vertebral body due to tumor or trauma.

Contraindications
Use of the coLigne ostaPek® VBR cages is contraindicated when there was recently or is an active infection, whether systemic or localized to the site of the proposed implantation. Any condition that precludes the possibility of fusion is relative contraindication. This includes but is not limited to: fever, mental illness, alcoholism or drug abuse, osteoporosis or osteopenia, obesity, pregnancy, foreign body sensitivity, extreme instability.
See also the warnings, precautions and possible adverse affects sections of this insert. The following are specific warnings, precautions, and adverse effects which should be understood by the surgeon and explained to the patient. General surgical risks should be explained to the patient prior to surgery.

Warnings
Potential risks identified with the use of the coLigne ostaPek® VBR cage, which may require additional surgery, include:
1. Implant fracture,
2. Loss of fixation,
3. Non-union,
4. Fracture of the vertebra,
5. Neurologic injury,
6. Vascular or visceral injury.
See also precautions and possible adverse events sections of the package insert for a complete list of potential risks.

Possible adverse events
Potential risks identified with the use of the coLigne ostaPek® VBR cages, which may require additional surgery, include but are not limited to:
1. Bending, fracture or loosening of implants,
2. Nonunion or delayed union,
3. Fracture of the vertebra,
4. Neurological, vascular or visceral injury,
5. Allergic reaction to a foreign body,
6. Infection,
7. Decrease in bone density due to stress shielding,
8. Pain, discomfort, or abnormal sensations due to the presence of the device,
9. Nerve damage due to surgical trauma,
10. Subsidence,
11. Bursitis,
12. Dural leak,
13. Paralysis,
14. Death.

Note
Additional surgery may be necessary to correct some of these anticipated adverse reactions.

Limited warranty
coLigne AG products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

Manufacturer
coLigne AG, Utoquai 43, CH-8008 Zurich
Phone +41 43 343 8000, Fax +41 43 343 8009
info@coligne.com, www.coligne.com

EU authorized representative
Axxos GmbH, Im Sägenlosh 3,
DE-78333 Stockach
eu_rep@axxos-medical.de

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Axxos GmbH, Im Sägentoh 3,